

# Instructions for Use

## Life Science Kits & Assays



### innuPREP Forensic Kit

**Order No.:**

845-KS-1051010 10 reactions

845-KS-1051050 50 reactions

845-KS-1051250 250 reactions

---

Publication No.: HB\_KS-1051\_e\_180328

---

This documentation describes the state at the time of publishing.  
It needs not necessarily agree with future versions. Subject to change!

Print-out and further use permitted with indication of source.

© Copyright 2018, Analytik Jena AG, AJ Innuscreen GmbH

---

**Manufacturer:**

AJ Innuscreen GmbH  
Robert-Rössle-Straße 10  
13125 Berlin  
Made in Germany!

---

**Distribution/Publisher:**

Analytik Jena AG  
Konrad-Zuse-Straße 1  
07745 Jena · Germany

Phone +49 3641 77 9400  
Fax +49 3641 77 767776  
[www.analytik-jena.com](http://www.analytik-jena.com)  
[info@analytik-jena.com](mailto:info@analytik-jena.com)

**CONTENTS**

|      |   |    |
|------|---|----|
| 1    | Introduction .....  | 2  |
| 1.1  | Intended use.....   | 2  |
| 1.2  | Notes on the use of this manual .....   | 3  |
| 2    | Safety precautions .....  | 4  |
| 3    | Storage conditions .....  | 5  |
| 4    | Functional testing and technical assistance .....   | 6  |
| 5    | Product use and warranty .....  | 6  |
| 6    | Kit components .....  | 7  |
| 7    | Product specifications.....   | 9  |
| 8    | GHS classification.....   | 10 |
| 8.1  | Hazard phrases .....  | 10 |
| 8.2  | Precaution phrases .....  | 11 |
| 8.3  | EU hazard statements .....  | 11 |
| 9    | Recommended steps before starting .....   | 12 |
| 10   | General procedure for DNA extraction.....   | 12 |
| 11   | Protocols .....   | 13 |
| 11.1 | Protocol 1: DNA isolation from buccal swab samples .....  | 13 |
| 11.2 | Protocol 2: DNA isolation from buccal swab samples from<br>different surfaces (cups, bottles, fingerprints etc.)..... | 16 |
| 11.3 | Protocol 3: DNA isolation from blood stains, saliva stains,<br>stamps and envelopes, etc.....                         | 19 |
| 11.4 | Protocol 4: DNA extraction from hair roots, barb hairs,<br>finger nails, etc.....                                     | 21 |
| 11.5 | Protocol 5: DNA extraction from cigarette butts .....   | 23 |
| 11.6 | Protocol 6: DNA extraction from chewing gum .....   | 25 |
| 11.7 | Protocol 7: DNA extraction from tissue samples .....  | 27 |
| 12   | Troubleshooting.....  | 29 |
| 13   | Related products .....  | 30 |









# 1 Introduction

## 1.1 Intended use

The kit has been designed as a tool for very fast and efficient isolation of genomic DNA from small amounts of different types of forensic samples like, hairs or hair roots; stains of blood, saliva or sperm; finger nails; cigarette butts; bubble gum; buccal swabs; stamps and envelopes as well as fingerprints on different surfaces. The kit is intended for use by professional users. The kit has been designed to be used for a wide range of different downstream application like amplification reactions and further analytical procedures. Diagnostic results generated using the extraction procedure in conjunction with diagnostic tests should be interpreted with regard to other clinical or laboratory results. To reduce irregularities in diagnostic results, internal controls for downstream applications should be used.

## 1.2 Notes on the use of this manual

For easy reference and orientation, the manual uses the following warning and information symbols as well as the shown methodology:

| Symbol  | Information   |
|---|---|
|    | <b>REF</b><br>Catalogue number.   |
|    | <b>Content</b><br>Contains sufficient reagents for <N> reactions.   |
|    | <b>Storage conditions</b><br>Store at room temperature or shown conditions respectively.  |
|    | <b>Consult instructions for use</b><br>This information must be observed to avoid improper use of the kit and the kit components. |
|  | <b>Expiry date</b>  |
|  | <b>Lot number</b><br>The number of the kit charge.  |
|  | <b>Manufactured by</b><br>Contact information of manufacturer.  |
|  | <b>For single use only</b><br>Do not use components for a second time.  |

The following systematic approach is introduced in the manual:

- The chapters and figures are numbered consecutively.
- A cross reference is indicated with an arrow (e.g. → „Notes on the use of this manual“ p. 3).
- Working steps are numbered.

## 2 Safety precautions

---

### NOTE

Read through this chapter carefully prior to guarantee your own safety and a trouble-free operation.

Follow all the safety instructions explained in the manual, as well as all messages and information, which are shown.

---

All due care and attention should be exercised in handling the materials and reagents contained in the kit. Always wear gloves while handling these reagents and avoid any skin contact! In case of contact, flush eyes or skin with a large amount of water immediately.



### FOR SINGLE USE ONLY!

This kit is made for single use only!

---

### ATTENTION!

Don't eat or drink components of the kit!

The kit shall only be handled by educated personnel in a laboratory environment!

---

If the buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries. This kit could be used with potential infectious samples. Therefore, all liquid waste must be considered as potentially infectious and must be handled and discarded according to local safety regulation.

Please observe the federal, state and local safety and environmental regulations. Follow the usual precautions for applications using extracted nucleic acids. All materials and reagents used for DNA or RNA isolation should be free of DNases or RNases.

---

**ATTENTION!**

Do not add bleach or acidic components to the waste after sample preparation!

---

**NOTE**

Emergency medical information in English and German can be obtained 24 hours a day from:

Poison Information Center, Freiburg / Germany

Phone: +49 (0)761 19 240.

---

For more information, please ask for the Safety Data Sheet (SDS).

### 3 Storage conditions

Store lyophilized **Proteinase K** at 4 °C to 8 °C! Divide dissolved **Proteinase K** into aliquots and store at -22 °C to -18 °C. Repeated freezing and thawing will reduce the activity dramatically!

All other components of the innuPREP Forensic Kit should be stored dry at room temperature (15 °C to 30 °C). When stored at room temperature, the kit is stable until the expiration date printed on the label on the kit box.

Before every use make sure that all components have room temperature. If there are any precipitates within the provided solutions dissolve these precipitates by careful warming.

For further information see chapter "Kit components" p.7.

## 4 Functional testing and technical assistance

The Analytik Jena AG guarantees the correct function of the kit for applications as described in the manual. This product has been produced and tested in an ISO 13485 certified facility.

We reserve the right to change or modify our products to enhance their performance and design. If you have any questions or problems regarding any aspects of the innuPREP Forensic Kit or other Analytik Jena AG products, please do not hesitate to contact us. For technical support or further information in Germany please dial +49 36 41 / 77 94 00. For other countries please contact your local distributor.

## 5 Product use and warranty

The kit is not designed for the usage of other starting materials or other amounts of starting materials than those, referred to in the manual (→ "Intended use" p. 2), (→ "Product specifications" p. 9). Since the performance characteristics of Analytik Jena AG kits have just been validated for the application described above, the user is responsible for the validation of the performance of Analytik Jena AG kits using other protocols than those described below. Analytik Jena AG kits may be used in clinical diagnostic laboratory systems after the laboratory has validated the complete diagnostic system as required by CLIA' 88 regulations in the U.S. or equivalents in other countries.

All products sold by Analytik Jena AG are subjected to extensive quality control procedures and are warranted to perform as described when used correctly. Any problems should be reported immediately.

---

### NOTE

The kit is for research use only!


---



## 6 Kit components

### IMPORTANT




Store lyophilized **Proteinase K** at 4 °C to 8 °C! Aliquot dissolved **Proteinase K** and store at -22 °C to -18 °C. Repeated freezing and thawing will reduce the activity dramatically!



### STORAGE CONDITIONS

All other components are stored at room temperature.

|                                | Σ<br>10                            | Σ<br>50                            | Σ<br>250                           |
|--------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <b>REF</b>                     | 845-KS-1051010                     | 845-KS-1051050                     | 845-KS-1051250                     |
| Lysis Solution TLS             | 5 ml                               | 25 ml                              | 120 ml                             |
| Binding Solution TBS           | 5 ml                               | 25 ml                              | 120 ml                             |
| Proteinase K                   | for 1 x 0.3 ml<br>working solution | for 1 x 1.5 ml<br>working solution | for 5 x 1.5 ml<br>working solution |
| Washing Solution HS<br>(conc.) | 3 ml                               | 15 ml                              | 70 ml                              |
| Washing Solution MS<br>(conc.) | 3 ml                               | 15 ml                              | 60 ml                              |
| Elution Buffer                 | 2 ml                               | 10 ml                              | 30 ml                              |
| Spin Filter                    | 10                                 | 50                                 | 5 x 50                             |
| Receiver Tubes                 | 40                                 | 4 x 50                             | 20 x 50                            |
| Elution Tubes                  | 10                                 | 50                                 | 5 x 50                             |
| Manual                         | 1                                  | 1                                  | 1                                  |

|                      |  10  |  50  |  250  |
|----------------------|---|---|--|
| <b>Initial steps</b> | <p><b>Proteinase K</b><br/>Dissolve Proteinase K by addition of 0.3 ml ddH<sub>2</sub>O, mix thoroughly and store as described above!</p> <p><b>Washing Solution HS</b><br/>Add 3 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> <p><b>Washing Solution MS</b><br/>Add 7 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> | <p><b>Proteinase K</b><br/>Dissolve Proteinase K by addition of 1.5 ml ddH<sub>2</sub>O, mix thoroughly and store as described above!</p> <p><b>Washing Solution HS</b><br/>Add 15 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> <p><b>Washing Solution MS</b><br/>Add 35 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> | <p><b>Proteinase K</b><br/>Dissolve Proteinase K by addition of 1.5 ml ddH<sub>2</sub>O, mix thoroughly and store as described above!</p> <p><b>Washing Solution HS</b><br/>Add 70 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> <p><b>Washing Solution MS</b><br/>Add 140 ml of 96-99.8 % ethanol to the bottle and mix thoroughly. Keep the bottle always firmly closed!</p> |

---

#### COMPONENTS NOT INCLUDED IN THE KIT

- RNase A (10 mg/ml); optional
- 1 M DTT solution
- 1.5 ml reaction tubes
- ddH<sub>2</sub>O for dissolving Proteinase K
- 96–99.8 % ethanol, non-denatured or methylated

## 7 Product specifications

### 1. Starting material:

- Swabs from different surfaces (e.g. cups, bottles, fingerprints)
- Blood samples
- Sperm samples
- Hair, hair roots or barb hairs
- Envelopes
- Finger nails
- Cigarette butts or paper
- Chewing gum

### 2. Time for isolation:

- Approximately 15 minutes after lysis step

---

#### NOTE

Using the kit for other kinds of forensic sample which are not described in the protocols, the selection of one of the described protocols is recommended.

---

---

#### NOTE

Optionally, for the isolation of DNA from forensic samples containing extremely low amounts of DNA it could be helpful to add **Carrier RNA** after lysis to the Binding Step.





We recommend to use Carrier RNA (e.g. Poly(A) RNA; Roche Diagnostics; No.108626). Dissolve the RNA in RNase-free water to obtain a solution of 1 µg/µl. Divide into aliquots and store at -22 °C to -18 °C.

Do not freeze and thaw the aliquots more than 3 times.

We recommend the addition of **1 µl Carrier RNA** per sample.

---

## 8 GHS classification

| Component                   | Hazard contents                                     | GHS Symbol   | Hazard phrases                | Precaution phrases  | EUH |
|-----------------------------|---|--|-------------------------------|---|-----|
| Binding Solution TBS        | Propan-2-ol<br>25-50 %                              |             | 225, 315,<br>318, 336,<br>411 | 101, 102, 103,<br>210,<br>303+361+353,<br>305+351+338,<br>310, 405, 501 |     |
|                             | Polyethylene glycol<br>octylphenol ether<br>25-50 % | <br>Danger  |                               |   |     |
| Washing Solution HS (conc.) | Guanidinium thiocyanate<br>50-100 %                 | <br>Danger  | 302, 314,<br>412              | 101, 102, 103,<br>260,<br>303+361+353,<br>305+351+338,<br>310, 405, 501 | 032 |
| Proteinase K                | Proteinase, engyodontium album<br>50-100 %          | <br>Danger | 315, 317,<br>319, 334,<br>335 | 101, 102, 103,<br>261, 280,<br>305+351+338,<br>342+311, 405,<br>501     |     |

### 8.1 Hazard phrases

|     |  |
|-----|--|
| 225 | Highly flammable liquid and vapour.  |
| 302 | Harmful if swallowed.  |
| 314 | Causes severe skin burns and eye damage.                                   |
| 315 | Causes skin irritation.  |
| 317 | May cause an allergic skin reaction.                                       |
| 318 | Causes serious eye damage.   |
| 319 | Causes serious eye irritation.   |
| 334 | May cause allergy or asthma symptoms or breathing difficulties if inhaled. |
| 335 | May cause respiratory irritation.  |
| 336 | May cause drowsiness or dizziness.   |
| 411 | Toxic to aquatic life with long-lasting effects.                           |

412 Harmful to aquatic life with long lasting effects.

## 8.2 Precaution phrases

- 101 If medical advice is needed, have product container or label at hand.
- 102 Keep out of reach of children.
- 103 Read label before use.
- 210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- 260 Do not breathe dust/fume/gas/mist/vapors/spray.
- 261 Avoid breathing dust/fume/gas/mist/vapors/spray.
- 280 Wear protective gloves/protective clothing/ eye protection/face protection.
- 310 Immediately call a POISON CENTER/doctor.
- 405 Store locked up.
- 501 Dispose of contents/container in accordance with local/regional/national/international regulations.
- 342+311 If experiencing respiratory symptoms: Call a POISON CENTER/doctor.
- 303+361+353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- 305+351+338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.


## 8.3 EU hazard statements

032 Contact with acids liberates very toxic gas.

## 9 Recommended steps before starting

- Ensure that the **Proteinase K / Washing Solution HS / Washing Solution MS** have been prepared according to the instruction (→ “Kit components” p. 7).
- Heat thermal mixer or water bath at 50 °C (optional 42 °C).
- Centrifugation steps should be carried out at room temperature.
- Avoid freezing and thawing of starting material.

## 10 General procedure for DNA extraction



## 11 Protocols

### 11.1 Protocol 1: DNA isolation from buccal swab samples

---

#### NOTE

To get maximum yield of DNA it is essential to leave the swab during the complete lysis time in the 1.5 ml tube. It is possible to cut the shaft of the swab, so that you can close the cap of the tube. The removal of the swab from the tube ahead of time will lead to a dramatically reduced final yield!

---

1. Place the swab into a 1.5 ml tube.

Add 400  $\mu$ l Lysis Solution TLS and 25  $\mu$ l Proteinase K, mix vigorously by pulsed vortexing for 5 seconds. Incubate at 50 °C for 10-15 minutes.

After lysis remove the swab from the tube and squeeze the swab on the wall of the tube to remove all Lysis Solution TLS from the swab.

---

#### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

#### NOTE

To remove RNA from the sample (if necessary) add 2  $\mu$ l of RNase A solution (10 mg/ml), vortex shortly and incubate for 5 minute at RT.

---

2. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.

---

### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

3. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

4. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
5. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
6. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
7. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **50–100 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step will increase the yield of extracted DNA.



**NOTE**

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

## 11.2 Protocol 2: DNA isolation from buccal swab samples from different surfaces (cups, bottles, fingerprints etc.)

---

### NOTE

To get maximum yield of DNA it is essential to leave the swab during the complete lysis time in the 1.5 ml tube. It is possible to cut the shaft of the swab, so that you can close the cap of the tube. The removal of the swab from the tube ahead of time will lead to a dramatically reduced final yield!

---

1. Place the swab into a 1.5 ml tube.  
Add **400 µl Lysis Solution TLS** and **25 µl Proteinase K**, mix vigorously by pulsed vortexing for 5 seconds. Incubate at 50 °C for 10-15 minutes.

After lysis remove the swab from the tube and squeeze the swab on the wall of the tube to remove all Lysis Solution TLS from the swab.

---

### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

### NOTE

To remove RNA from the sample (if necessary) add 2 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 5 minute at RT.

---

2. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.

---

**NOTE**

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

3. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

**NOTE**

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

4. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
5. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
6. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
7. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **50–100 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step will increase the yield of extracted DNA.

---

### **NOTE**

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

### 11.3 Protocol 3: DNA isolation from blood stains, saliva stains, stamps and envelopes, etc.

1. Cut the material containing the stains into small pieces and transfer it into a 1.5 ml reaction tube.

Add **400 µl Lysis Solution TLS** and **25 µl Proteinase K**.

For semen stains add **30 µl 1 M DTT** (not provided) to the Lysis Solution/ Proteinase K mix.

Mix vigorously by pulsed vortexing for 5 seconds. Incubate at 50 °C for at least 2 hours.

---

#### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

2. Centrifuge the 1.5 ml tube at 11,000 x g (~11,000 rpm) for 1 minute to spin down unlysed material.

Transfer the supernatant into another 1.5 ml tube.

3. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.
- 

#### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

4. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

5. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
6. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
7. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
8. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **30 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step will increase the yield of extracted DNA.

---

### NOTE

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

#### 11.4 Protocol 4: DNA extraction from hair roots, barb hairs, finger nails, etc.

1. Cut the material into small pieces and transfer it into a 1.5 ml reaction tube.  
Add 400 µl Lysis Solution TLS and 25 µl Proteinase K.
2. Add 30 µl 1 M DTT (not provided) to the Lysis Solution / Proteinase K mix. Mix vigorously by pulsed vortexing for 5 seconds.  
Incubate at 50 °C for at least 2 hours (incubation overnight at 42 °C is also possible).

**Assure that the hair roots are in the Lysis Solution during the lysis time!**

---

#### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

3. Centrifuge the 1.5 ml tube at 11,000 x g (~11,000 rpm) for 1 minute to spin down unlysed material.  
Transfer the supernatant into another 1.5 ml tube.
  4. Add 400 µl Binding Solution TBS to the lysed sample, mix by brief vortexing or by pipetting up and down several times.
- 

#### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

5. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

6. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
7. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
8. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
9. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **30 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step will increase the yield of extracted DNA.

---

### NOTE

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---



### 11.5 Protocol 5: DNA extraction from cigarette butts

1. Remove of a small piece (3–5 mm) of the brown filter paper or of a part of the filter and place the material in a 1.5 ml reaction tube. Add **400 µl of Lysis Solution TLS** and **25 µl of Proteinase K**, mix vigorously by pulsed vortexing for 5 seconds.

Incubate at 50 °C for at least 2 hours (incubation overnight at 42 °C is also possible).

---

#### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

2. Centrifuge the 1.5 ml tube at 11,000 x g (~11,000 rpm) for 1 minute to spin down unlysed material.  
Transfer the supernatant into another 1.5 ml tube.
3. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.

---

#### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

4. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
- 

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

5. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
  6. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
  7. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
  8. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **30 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step (with other 30 µl Elution Buffer) will increase the yield of extracted DNA.
- 

### NOTE

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

## 11.6 Protocol 6: DNA extraction from chewing gum

1. Cut a part of the chewing gum into small pieces and place the material into a 1.5 ml reaction tube.  
Add **400 µl Lysis Solution TLS** and **25 µl Proteinase K**, mix vigorously by pulsed vortexing for 5 seconds.

Incubate at 50 °C for at least 2 hours (incubation overnight at 42 °C is also possible).

---

### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

2. Centrifuge the 1.5 ml tube at 11,000 x g (~11,000 rpm) for 1 minute to spin down unlysed material.  
Transfer the supernatant into another 1.5 ml tube.
3. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.

---

### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

4. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

5. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
6. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
7. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
8. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **30 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step (with other 30 µl Elution Buffer) could increase the yield of extracted DNA.

---

### NOTE

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

### 11.7 Protocol 7: DNA extraction from tissue samples

1. Cut the material (fresh or frozen) into small pieces and transfer it into a 1.5 ml reaction tube.

Add **400 µl Lysis Solution TLS** and **25 µl Proteinase K**, mix vigorously by pulsed vortexing for 5 seconds.

Incubate at 50 °C for at least 1 hour (incubation overnight at 42 °C is also possible).

---

#### NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

---

2. Centrifuge the 1.5 ml tube at 11,000 x g (~11,000 rpm) for 1 minute to spin down unlysed material.

Transfer the supernatant into another 1.5 ml tube.

---

#### NOTE

To remove RNA from the sample (if necessary) add 2 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 5 minute at RT.

---

3. Add **400 µl Binding Solution TBS** to the lysed sample, mix by brief vortexing or by pipetting up and down several times.
- 

#### NOTE

It is important that the sample and the Binding Solution TBS are mixed vigorously to get a homogeneous solution.

---

4. Apply the sample to the Spin Filter located in a Receiver Tube. Close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.

---

### NOTE

If the solution has not completely passed through the Spin Filter, centrifuge again at higher speed or prolong the centrifugation time.

---

5. Open the Spin Filter and add **500 µl Washing Solution HS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
6. Open the Spin Filter and add **750 µl Washing Solution MS**, close the cap and centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. Discard the Receiver Tube with the filtrate. Place the Spin Filter into a new Receiver Tube.
7. Centrifuge at 11,000 x g (~11,000 rpm) for 2 minutes to remove all traces of ethanol. Discard the Receiver Tube.
8. Place the Spin Filter into an Elution Tube. Carefully open the cap of the Spin Filter and add **30 µl Elution Buffer**. Incubate at room temperature for 1 minute. Centrifuge at 11,000 x g (~11,000 rpm) for 1 minute. A second elution step (with other 30 µl Elution Buffer) could increase the yield of extracted DNA.

---

### NOTE

The DNA can be eluted with a lower or a higher volume of Elution Buffer (depends on the expected yield of genomic DNA). Elution with lower volumes of Elution Buffer increases the final concentration of DNA. Store the extracted DNA at 4 °C to 8 °C. For long time storage placing at -18 °C to -22 °C is recommended.

---

## 12 Troubleshooting

| Problem / probable cause                             | Comments and suggestions   |
|--|--|
| <b>Clogged Spin Filter</b>                           |  |
| Insufficient lysis and/or too much starting material | Increase lysis time.<br>Increase centrifugation speed.<br>After lysis centrifuge the lysate to pellet un-lysed material.<br>Reduce amount of starting material.  |
| <b>Low amount of extracted DNA</b>                   |  |
| Insufficient lysis                                   | Increase lysis time!<br>Reduce amount of starting material. Overloading reduces yield!   |
| Incomplete elution                                   | Prolong the incubation time with Elution Buffer to 5 minutes or repeat elution step once again.<br>Take a higher volume of Elution Buffer.   |
| Insufficient mixing with Binding Solution TBS        | Mix sample with Binding Solution TBS by pipetting or by vortexing prior to transfer of the sample onto the Spin Filter.  |
| <b>Low concentration of extracted DNA</b>            |  |
| Too much Elution Buffer was used in the elution step | Elute the DNA with lower volume of Elution Buffer  |
| <b>Degraded or sheared DNA</b>                       |  |
| Incorrect storage of starting material               | Ensure that the starting material is frozen immediately after taking in liquid nitrogen or at -18 °C to -22° C! For long time storage continuously store at -78 °C to -82° C! Avoid thawing of the material. |
| Old starting material                                | Old material often contains degraded DNA. Repeat with fresh material.  |
| <b>RNA contamination</b>                             |  |
| Extracted DNA is contaminated with RNA               | Perform an RNase A digestion.  |

## 13 Related products

| Name  | Amount     | Order No.      |
|---|------------|----------------|
| <b>Products for PCR &amp; Gel Electrophoresis</b> |            |                |
| innuPREP DOUBLEpure Kit                           | 10 rxn     | 845-KS-5050010 |
|   | 50 rxn     | 845-KS-5050050 |
|   | 250 rxn    | 845-KS-5050250 |
| innuPREP Gel Extraction Kit                       | 10 rxn     | 845-KS-5030010 |
|   | 50 rxn     | 845-KS-5030050 |
|   | 250 rxn    | 845-KS-5030250 |
| innuPREP PCRpure Kit                              | 10 rxn     | 845-KS-5010010 |
|   | 50 rxn     | 845-KS-5010050 |
|   | 250 rxn    | 845-KS-5010250 |
| innuTaq DNA Polymerase (5 U/μl)                   | 500 U      | 845-EZ-1000500 |
| 50x inNucleotide Mix (1.5 mM)                     | 2x 0.5 ml  | 845-AS-9000100 |
| inNucleotide Set (100 mM)                         | 4x 0.25 ml | 845-AS-1100250 |
| innuMIX rapidPCR MasterMix                        | 100 rxn    | 845-AS-1600100 |
|   | 200 rxn    | 845-AS-1600200 |
| innuMIX Standard PCR MasterMix                    | 100 rxn    | 845-AS-1700100 |
|   | 200 rxn    | 845-AS-1700200 |
| innuMIX Green PCR MasterMix                       | 100 rxn    | 845-AS-1400100 |
|   | 200 rxn    | 845-AS-1400200 |
| innuSTAR 100 bp DNA Ladder Express                | 500 μl     | 845-ST-1010100 |
|   | 5x 500 μl  | 845-ST-1010500 |
| innuSTAR 1 kb DNA Ladder Express                  | 500 μl     | 845-ST-1020100 |
|   | 5x 500 μl  | 845-ST-1020500 |
| 6x Loading Dye Bromophenol Blue                   | 3x 1.0 ml  | 845-ST-3010003 |
|   | 6x 1.0 ml  | 845-ST-3010006 |





#### Headquarters

---

Analytik Jena AG  
Konrad-Zuse-Str. 1  
07745 Jena · Germany

Phone +49 3641 77 70  
Fax +49 3641 77 9279  
info@analytik-jena.com  
www.analytik-jena.com

Pictures: Analytik Jena AG  
Subject to changes in design and scope of delivery as well as further technical development!