

Instructions for Use

Life Science Kits & Assays







innuTaq HOT-A DNA Polymerase

1 Product specifications

The innuTaq HOT-A DNA Polymerase provides improved specificity and sensitivity when amplifying low-copy-number targets in complex backgrounds. The polymerase activity is blocked at ambient temperature with a polymerase specific aptamere and switched on at the onset of the initial denaturation. The thermal activation prevents the extension of non specifically annealed primers and primer-dimer formation at low temperatures during PCR set-up.

It catalyses the polymerisation of nucleotides into duplex DNA in 5' → 3' direction in the presence of magnesium. It also possesses a 5' → 3' polymerisation – dependent exonuclease replacement activity but lacks a 3' → 5' exonuclease activity

Delivered components

Components	Description	Amount
 innuTaq HOT-A DNA Polymerase	Concentration: 5 U/μl Enzyme storage buffer: 20 mM Tris-HCl, 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5 % Tween 20, 0.5 % Nonidet P40, 50 % Glycerol, pH 8.0. I	500 U
 10x Hot Start Buffer complete	200 mM Tris-HCl (pH 8.5), 500 mM KCl, 15 mM MgCl ₂	1500 μl
 10x Hot Start Buffer without MgCl ₂	200 mM Tris-HCl (pH 8.5), 500 mM KCl	1500 μl
 Mg ²⁺ Solution	25 mM MgCl ₂	1500 μl

2 Product and order number

Name	Amount	Order-no.
innuTaq HOT-A DNA Polymerase	500 U	845-EZ-3000500

3 Storage conditions

Store innuTaq HOT-A DNA Polymerase at -22 to -18 °C in a freezer with constant temperature conditions.

When stored as recommended, the innuTaq HOT-A DNA Polymerase is stable until the expiration date printed on the label on the kit box.

4 Quality data and unit definition

Activity and stability tested by low copy PCR, human DNA contamination and activity of DNase and RNase are not detected. Polymerization activity at 25 °C is not detected.

One unit of enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides (dNTP's) into a polynucleotide fraction in 30 minutes at 70 °C.

5 Safety precautions

The assay shall only be handled by educated personal in a laboratory environment. The compliance with the specified procedure is absolutely mandatory when performing this assay.

Reagents should be stored in their original containers at the indicated temperatures. Do not replace individual components with those from different batches or test assays. Note the indicated expiration dates.

Do not eat, drink or smoke while performing the assay.

Wear protective clothing and safety gloves.

All samples and test materials should be handled and disposed of as infectious material, in accordance with regulatory requirements.

Reagent containers that have not come in contact with potentially infectious material may be disposed of along with ordinary laboratory waste.

Store the reagents used for performing PCR separately from DNA templates and amplification products.

6 Reagent preparation

- After thawing gently vortex and briefly centrifuge all solutions.
- Prepare a mix of following components.

Reagent	Volume (1 rxn)
10x Hot Start Buffer (complete or without MgCl ₂)	5 µl
25 mM MgCl ₂ Solution (if required)	3-5 µl
12.5 mM dNTP Mix	1 µl
Forward Primer	0.2 - 1 µM
Reverse Primer	0.2 - 1 µM
innuTaq HOT-A DNA Polymerase (5 U/µl)	0.2 - 0.5 µl
Template DNA	1 - 100 ng/µl (≤ 1 µg)
PCR-grade H ₂ O	add to a final vol. of 50 µl
Total volume	50 µl

7 PCR conditions

Step	Cycles	Profile	Temperature	Retention time
1	1	Initial denaturation	95 °C	120 s
		Denaturation	95 °C	30 - 60 sec
2	25-40	Annealing	50 - 68 °C	30 - 60 sec
		Elongation	72 °C	1 - 4 min
3	1	Final elongation	72 °C	5 - 10 min

Note: Annealing temperature should be 2 - 6 °C lower than melting temperature of primer.

8 Hints and Notes

- Gently vortex and briefly centrifuge all solutions after thawing
- After pipetting mix the components of the reaction mix by gently vortexing and briefly centrifuge for a few seconds to collect the mixture at the bottom of the tube.
- The enzyme allows the PCR set up at ambient temperature without non-specific annealing and extension.

Reaction conditions (incubation temperatures and times, concentrations of template DNA, primers, magnesium ions and enzyme) are depending on the used template and primers.

The optimal Mg^{2+} concentrations vary between 1 - 4 mM and have to be determined empirically. However, most applications work at the standard concentration of 1.5 mM Mg^{2+} . Advanced applications on genomic DNA require higher Mg^{2+} concentrations (2 - 3 mM) adjustable by using the separate 25 mM $MgCl_2$ Solution supplied with the set.

9 Related products

Product	Order Number
50x inNucleotide Mix	845-AS-9000100
inNucleotide Set (100 mM)	845-AS-1100250
6x Loading Dye Orange G	845-ST-4010006
6x Loading Dye Bromophenol Blue	845-ST-3010006
innuSTAR 100 bp DNA Ladder Express	845-ST-1010500
innuSTAR 1kb DNA Ladder Express	845-ST-1020500

Publication No.: HB_EZ-3000_e_170817

This documentation describes the state at the time of publishing. It needs not necessarily agree with future versions. Subject to change!

Expression and further use permitted with indication of source. © 2017 Analytik Jena AG,
AJ Innuscreen GmbH

Manufacturer:

AJ Innuscreen GmbH
Robert-Rössle-Strasse 10
13125 Berlin · Germany

Distribution/Publisher:

Analytik Jena AG
Konrad-Zuse-Strasse 1
07745 Jena · Germany

Telefon +49 3641 77 70
Telefax +49 3641 77 9279
info@analytik-jena.com
www.analytik-jena.com