








## innuPREP DNA / RNA Mini Kit

### Protocol 1: DNA and RNA extraction from tissue samples


Recommended steps before starting

- Prepare Washing Solution HS and Washing Solution LS according to the instruction

1. Starting material	<ul style="list-style-type: none"> <li>Tissue</li> </ul>	<ul style="list-style-type: none"> <li>Max. 20 mg</li> </ul>	
2. Homogenization and lysis	<ul style="list-style-type: none"> <li>Homogenizer e.g. SpeedMill</li> </ul> <p>or</p> <ul style="list-style-type: none"> <li>Liquid nitrogen</li> </ul>	    	<ul style="list-style-type: none"> <li>Add frozen starting material to homogenizer tube</li> <li>Add 450 µl Lysis Solution RL and homogenize</li> <li>Add sample to a 1.5 ml tube</li> <li>Grind starting material to fine powder under liquid nitrogen</li> <li>Add sample to a 1.5 ml tube</li> <li>Add 450 µl Lysis Solution RL and lyse sample under continuous shaking</li> <li>Centrifuge: max. speed; 1 min</li> </ul>
3. Binding of DNA New Receiver Tube			<ul style="list-style-type: none"> <li>Spin Filter D to Receiver Tube</li> <li>Add supernatant to Spin Filter D</li> <li>10,000 x g (~12,000 rpm): 2 min</li> <li>Spin Filter D to Receiver Tube</li> <li><b>Don't</b> discard the filtrate!</li> </ul>
4. Binding of RNA New Receiver Tube			<ul style="list-style-type: none"> <li>Add equal volume 70 % ethanol (approx. 400 µl) to filtrate</li> <li>Spin Filter R to Receiver Tube</li> <li>Add filtrate to Spin Filter R</li> <li>10,000 x g (~12,000 rpm): 2 min</li> <li>Spin Filter R to Receiver Tube</li> </ul>

## 7. Washing of Spin Filter D and R


Re-use Receiver Tube



- Add 500 µl Washing Solution HS
- 10,000 x g (~12,000 rpm): 1 min
- Add 700 µl Washing Solution LS
- 10,000 x g (~12,000 rpm): 1 min

## 8. Remove Ethanol of Spin Filter D and R

Re-use Receiver Tube



- Discard filtrate
- Spin Filter D to Receiver Tube
- Spin Filter R to Receiver Tube
- Centrifuge: max speed, 2 min

## 9. Elution of Spin Filter D and R



- Add Spin Filter to an Elution Tube
- Add 100–200 µl Elution Buffer
- Incubation: 1 min @ RT
- 11,000 x g (~11,000 rpm): 1 min

### Order No.:

845-KS-2080010

10 reactions

845-KS-2080050

50 reactions

845-KS-2080250

250 reactions

This documentation describes the state at the time of publishing. It needs not necessarily agree with future versions. Subject to change!

*Expression and further use permitted with indication of source. © 2017 Analytik Jena AG, AJ Innuscreen GmbH*

Analytik Jena AG  
Konrad-Zuse-Straße 1  
07745 Jena · Germany

Phone +49 3641 77 9400  
Fax +49 3641 77 76 7776  
info@analytik-jena.com  
www.analytik-jena.com

Manufacturer:  
AJ Innuscreen GmbH  
Robert-Rössle-Straße 10  
13125 Berlin · Germany