









## innuPREP Plasmid Mini Kit

### Protocol 1: Isolation of plasmid DNA from 0.5–5.0 ml bacterial culture

Recommended steps     ■ Prepare Washing Solution B according to the instruction before starting

1. Starting material	<i>E. coli</i> overnight culture	■ 0.5–5.0 ml
2. Pellet cells		<ul style="list-style-type: none"> <li>■ Centrifuge: max. speed, 1 min</li> <li>■ Discard supernatant completely</li> </ul>
3. Resuspend cells		<ul style="list-style-type: none"> <li>■ Add 250 µl Resuspension Buffer</li> <li>■ Vortex</li> </ul>
4. Lysis Don't vortex!		<ul style="list-style-type: none"> <li>■ Add 250 µl Lysis Buffer</li> <li>■ Mix: invert tube 6–8 times</li> <li>■ Lysis time: ≤ 5 min</li> </ul>
5. Neutralization	 	<ul style="list-style-type: none"> <li>■ Add 350 µl Neutralization Buffer</li> <li>■ Mix: invert tube 6–8 times</li> <li>■ Centrifuge: max speed, 8 min</li> </ul>
6. Binding of DNA	 	<ul style="list-style-type: none"> <li>■ Add Spin Filter to Receiver Tube</li> <li>■ Add clarified sample to Spin Filter</li> <li>■ 11,000 x g (~11,000 rpm): 1 min</li> </ul>
7. Washing Re-use Receiver Tube	 	<ul style="list-style-type: none"> <li>■ Add 500 µl Washing Solution A</li> <li>■ 11,000 x g (~11,000 rpm): 1 min</li> <li>■ Add 700 µl Washing Solution B</li> <li>■ 11,000 x g (~11,000 rpm): 1 min</li> </ul>

8. Remove  
Ethanol

Re-use Receiver Tube



- Discard filtrate
- Add Spin Filter to Receiver Tube
- Centrifuge: max speed, 2 min

9. Elution



- Add Spin Filter to a 1.5 ml reaction tube
- Add 50–100 µl Elution Buffer P
- Incubation: 1 min @ RT
- 11,000 x g (~11,000 rpm): 1 min

<b>Order No.:</b>	Spin Filter without cap
	845-KS-5040010      10 reactions
	845-KS-5040050      50 reactions
	845-KS-5040250      250 reactions
	845-KS-5040500      500 reactions
	Spin Filter with cap
	845-KS-5041010      10 reactions
	845-KS-5041050      50 reactions
	845-KS-5041250      250 reactions
	845-KS-5041500      500 reactions

This documentation describes the state at the time of publishing.  
It needs not necessarily agree with future versions. Subject to change!

*Expression and further use permitted with indication of source. © 2017 Analytik Jena AG, AJ Innuscreen GmbH*

Analytik Jena AG  
Konrad-Zuse-Strasse 1  
07745 Jena · Germany

Phone +49 3641 77 70  
Fax +49 3641 77 9279  
info@analytik-jena.com  
www.analytik-jena.com

**Manufacturer:**  
AJ Innuscreen GmbH  
Robert-Rössle-Strasse 10  
13125 Berlin · Germany