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
Table of Contents

1	Basic information.....	1
1.1	User manual notes	1
	Conventions.....	1
1.2	Intended use	1
2	Safety instructions	2
3	Technical data	3
4	Design of the Double Furnace Solids Module	4
	Components for sample preparation	4
	Combustion system	4
	Accessories	4
5	Installing the Double Furnace Solids Module	5
	Filling the combustion tube of the solids module	5
	Installing the solids module	6
6	Operation of the Double Furnace Solids Module	8
	Preparing for the measurement.....	8
	Tempering the sample boats.....	8
	Performing the analysis.....	8
7	Maintenance and care for the Double Furnace Solids Module.....	10
8	Removing the Double Furnace Solids Module	11

Index of Figures

Fig. 1	Double Furnace Solids Module	1
Fig. 2	Design of the Double Furnace Solids Module	4
Fig. 3	Inserting a sample boat in the Double Furnace Solids Module	9

1 Basic information

1.1 User manual notes

The Double Furnace Solids Module is a special reactor and a lock with manual feed for the analysis of small quantities of solid samples and can extend the combustion systems of the analyzers multi N/C 2100S / 2100. Digestion temperatures of up to 950 °C are achieved. The digestion is supported by catalysts.


Fig. 1 Double Furnace Solids Module

CAUTION

Carefully read the above-mentioned documents before connecting the Double Furnace Solids Module to the multi N/C 2100/2100S and using it. Observe in particular the notes in the chapters "Safety notes" of the user manuals.

Conventions

Instructions for action which occur in chronological order are numbered and combined into action units and furnished with the corresponding results.

Lists which are not in chronological order are shown as itemized lists, sub-listings as bullet points.

Safety notes are indicated by pictographs and signal words. The type and source of the danger as well as the consequences are stated together with notes on preventing the danger. The meaning of the pictographs and signal words used is explained in the chapter "Safety notes".

The elements of the control and analysis software are indicated as follows:

- Menu commands, buttons, options etc. are formatted in bold.
- Menu commands of a command sequence are separated by slashes (/), e.g. **Method / Edit**.
- Window names are shown in bold italics, e.g. ***Measurement*** window.
- Buttons are marked with square brackets, e.g. **[Save]**.

1.2 Intended use

The Double Furnace Solids Module only must be used for the methods for supplying solids to the combustion system of the multi N/C 2100/2100S as described in this operating manual. Any other use is not as intended!

2 Safety instructions

IMPORTANT

This user manual is only valid in conjunction with the following documents:

- User manual for multi NC 2100/2100S
- Software manual of the control and analysis software multiWin

Observe in particular the notes in the chapters "Safety notes" of these operating manuals. The notes included there apply correspondingly also to the Double Furnace Solids Module without any limitations.

Special dangers implied in the Double Furnace Solids Module will be indicated separately in the single chapters.

Observe all safety notes listed in this operating manual and observe all messages and notes displayed by the control and analysis software multiWin on the monitor.

Please observe the following notes when connecting and operating the Double Furnace Solids Module:

CAUTION

There is a risk of burning on the combustion furnace! Only remove the combustion tube when the device is cold or allow the device to cool down sufficiently!

Before switching off set the furnace temperature in multiWin to 20 °C and exit multiWin. Otherwise there is a risk of burns when checking the system for leaks after installation!

3 Technical data

Digestion temperature	up to 950 °C
Catalyst	CeO ₂ (special catalyst for multi N/C)
Sample volume	max. 500 mg (dependent on the sample material)
Sample Feed	Boats via manual lock
Carrier gas supply	Oxygen (min. 4.5) preliminary pressure 4 – 6 bar

4 Design of the Double Furnace Solids Module

The Double Furnace Solids Module consists of the following main components:

- Components for sample preparation
- Combustion system
- Accessories

The solids module is connected via an adapter at the mount of the furnace lock to the analyzer multi N/C 2100/2100S.


Fig. 2 Design of the Double Furnace Solids Module

Components for sample preparation

The furnace lock is mounted on the side opening of the combustion tube. For the solids module the samples are supplied via sample boats. The sample amount is max. 500 mg (dependent on the sample material). The sample boat is inserted into the solids module via the furnace lock. It is opened and closed manually via an interlock.

Combustion system

The use of the Double Furnace Solids Module in the analyzer multi N/C 2100/2100S requires the use of a combustion furnace for vertical and horizontal operation. The combined combustion furnace permits the operation with a vertically installed combustion tube or the operation with a horizontally installed solids module.

The horizontally installed combustion tube of the solids module is made from quartz glass. The furnace lock with manual feed is mounted to the side opening of the combustion tube. The gas transfer hose is connected to the combustion tube adapter and linked to the condensation coil in the analyzer via a fork clamp.

The double-walled combustion tube is filled with catalyst and auxiliary material. As catalyst in the Double Furnace Solids Module the special catalyst for multi N/C (CeO_2) with a reaction temperature of up to 950 °C (standard: 900 °C) is used as default.

Accessories

The following accessories are required to operate the Double Furnace Solids Module:

- Connection hoses
- Tools

5 Installing the Double Furnace Solids Module

The Double Furnace Solids Module can only be installed if a combined combustion furnace for horizontal and vertical operation has been installed.

Filling the combustion tube of the solids module

CAUTION!

Alkaline salts (hand perspiration) causes crystallization in the quartz glass when heating combustion tubes in the combustion furnace. This shortens the service life of the combustion tube.

Where possible do not touch the cleaned combustion tube with your hands for filling. Wear protective gloves during the filling of the combustion tube.

Only fill fully dried combustion tubes. If necessary, dry the combustion tube before filling.




Wipe any finger marks on the combustion tube with a cloth wetted with pure alcohol.

	<ol style="list-style-type: none"> 1. To fill rotate the adapter (4) of the combustion tube upwards. 2. Insert the quartz wool into the combustion tube, carefully push it down with a glass rod and compact it (2). 3. Carefully fill approx. 60 g CeO_2 (special catalyst multi N/C) via the adapter into the sleeve of the combustion tube (1). 4. Seal the adapter of the combustion tube with some quartz wool (3). <p><u>Note:</u></p> <p>The quartz wool is used to retain the catalyst. Seal the adapter to make sure that no catalyst can enter the gas path. However, do not compact the quartz wool too densely.</p>
	<ol style="list-style-type: none"> 5. Open the furnace lock: <ul style="list-style-type: none"> – Pull the interlock (1) upwards. – Pull the bracket (2) out of its anchorage. – Pull the furnace lock open.

	<ol style="list-style-type: none"> 6. Release the three hexagon socket screws with the corresponding angular screwdriver by half a turn. Do not fully unscrew the hexagon socket screws. 7. Slide the filled combustion tube up to the stop at the inner ring into the Double Furnace Solids Module. The adapter of the gas outlet must then point down. 8. Tighten the hexagon socket screws tightly. 9. Close the furnace lock.
--	---

Installing the solids module

<ol style="list-style-type: none"> 1. Remove the combustion tube for vertical operation (see section "Removing the combustion tube" in the operating manuals of the analyzers). 	
	<ol style="list-style-type: none"> 2. Remove the sealing plug from the horizontal opening of the combustion furnace and place it onto the vertical opening of the combustion furnace.
	<ol style="list-style-type: none"> 3. Attach the holding plate with the four knurled head screws to the angular profiles in front of the horizontal opening of the combustion furnace.
	<ol style="list-style-type: none"> 4. Insert the solids module into the horizontal opening of the combustion furnace. The adapter of the combustion tube points down. 5. Attach the solids module with three hexagon socket screws to the holding plate.
	<ol style="list-style-type: none"> 6. Attach the gas transfer hose to the adapter of the combustion tube.

	<ol style="list-style-type: none"> 7. Combine the gas transfer hose and the inlet of the condensation coil (spherical joint). 8. Secure the spherical joint with the fork clamp and tighten the knurled head screw of the gas clamp finger-tight.
	<ol style="list-style-type: none"> 9. Attach the gas inlet hose to the carrier gas outlet of the gas box by way of FAST connectors. 10. Screw the other end of the hose to the end of the Double Furnace Solids Module using the Fingertight connection.
	<ol style="list-style-type: none"> 11. Connect the plug-in connector of the analyzer to the left-hand side of the Double Furnace Solids Module. 12. Attach the protective conductor connection and close the side wall.

6 Operation of the Double Furnace Solids Module

Preparing for the measurement

IMPORTANT

Oxygen (4.5) must be used as carrier gas when using the Double Furnace Solids Module.

Switch on the analyzer as follows:

1. Before switching on the analyzer check that the solids module has been tightly installed and connected as specified.
2. Switch on the analyzer.
The carrier gas flow is automatically adjusted to 400 ml/min.
3. Carry out the system leak test (see section "Checking the system for tightness" in the operating manual of the analyzer).
4. Heat up the combustion furnace. To do so set the furnace temperature in the window **Options** / tab **Device components** to the desired temperature (e.g. 900 °C).

Tempering the sample boats

IMPORTANT

The sample boats might be contaminated. Temper the sample boats prior to analyzing standards and samples. Tempering is performed by an "empty measurement".

After tempering do not touch the sample boat by hand any more. Use e.g. a clean petri dish to store the sample boats and a clean pair of tweezers to transport them.

Performing the analysis

Note:

When using the solids module only measurements with manual sample feed are possible.

Prepare the analysis and start the measurement as follows:

1. Using the menu command **Method / New** create a new method or load an existing method.

To do so open via the menu command **Method / Load** the database window **Method selection**, highlight the desired method and confirm the selection by clicking the button **[OK]**.

Note:
When creating a new method the option **Horizontal furnace** must be selected.
2. In the window **System state** select the manual sample supply by clicking the button **[manual]**.
✓ This is followed by the initialization of the analyzer.
3. In the window **System state** check the following entries:
 - Visual bank - OK
 - Gas flow - OK
 - Temperature - OK

Note:

If one of the entries is incorrect (shown in red), carry out a fault finding exercise in accordance with the notes in the section "Fault removal" in the operating manual of the analyzer.

4. In the window **Options** / tab **FreeStrings** you may store information about the sample.
5. Click on **[Start measurement]** or open the menu command **Measurement / Start Measurement**.
The window **Measurement start** opens.
6. Enter the **Sample ID** and, if applicable, a name for the analysis table.
7. Define the **Sample type** for the solid sample to be measured.
8. If applicable, enter explanations about the measurement via the button **[Comment]**.
9. With **[Start ►]** open the window **Measurement**.
10. Start the measurement by clicking the button **[Start F2]** and follow the instructions of the control and analysis software.
11. In the window **Sample** enter the sample volume and then exit with **[OK]**.
12. After a prompt by the information window insert the sample (in the sample boat) into the furnace lock.
 - Open the furnace lock.
 - Insert the sample boats into the furnace lock and hang the eye of the sample boat into the hook of the feeder.

Note:

The lock should not be closed, otherwise highly volatile substances are not detected.


Fig. 3 Inserting a sample boat in the Double Furnace Solids Module

13. Exit the information window by clicking on the button **[OK]**.
14. Follow the prompts of another information window:
 - Close the lock.
 - Slide the sample to the combustion furnace.

Note:

If multiple detection has been agreed in the method, the process restarts at item (10).

- ✓ **At the end of the measurement the results appear in the analysis report or in the selected analysis table.**

7 Maintenance and care for the Double Furnace Solids Module

Combustion tube	
Maintenance task	Maintenance interval
Check for cracks and damages	3 months
Check catalyst and replace if necessary	as required, after corresponding message in multiWin at the latest
Clean combustion tube	as required, but no later than after 12 months

Furnace lock	
Maintenance task	Maintenance interval
Check for leaks	3 months
Replace the sealing ring of the lock	as required

8 Removing the Double Furnace Solids Module


CAUTION


There is a risk of burning on the combustion furnace! Only remove the combustion tube when the device is cold or allow the device to cool down sufficiently!

Before switching off set the furnace temperature in multiWin to 20 °C and exit multiWin. Otherwise there is a risk of burns when checking the system for leaks after installation!


1. In the window *Options* / tab *Device components* set a **furnace temperature** of 20 °C.
2. Exit the software multiWin.
3. Switch the analyzer off on the main switch, pull the mains plug from the mains outlet and disconnect the gas supply.



4. Disconnect the plug-in connector on the left-hand side of the solids module.

5. Pull the gas inlet hose out of the FAST connector at the carrier gas outlet of the gas box and remove the FAST connector from the terminal.
6. Unscrew the other end of the hose from the solids module.

7. Remove the fork clamp at the spherical joint between the gas transfer hose and the condensation coil inlet.

	<p>8. Loosen the four knurled head screws at the holding plate and pull the solids module out of the opening of the combustion furnace.</p> <p><i>Note:</i> The gas transfer hose and holding plate may remain connected to the solids module. This makes it easier to install the solids module the next time.</p> <p><i>Caution:</i> Do not unscrew the angle profiles from the furnace! They are preadjusted and guarantee the correct position of the solids module during installation.</p>
	<p>9. Remove the sealing plug from the vertical opening of the combustion furnace and place it onto the horizontal opening of the combustion furnace.</p>
<p>10. Install the combustion tube for vertical operation (see section "Removing the combustion tube" in the operating manual of the analyzer).</p>	

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Trademarks multi N/C and multiWin are registered trademarks of Analytik Jena AG in Germany.
Microsoft and Windows are registered trademarks of Microsoft Corp.
© Copyright 2011, Analytik Jena AG