

We Change the Way to Prep SmartExtraction



analytikjena
An Endress+Hauser Company

SmartExtraction Kicks off a New Era

Developed by Analytik Jena, the technology is the first of its kind in the world to facilitate simple automation with conventional liquid-handling systems², as well as requiring minimal equipment when running manual processes.

SmartExtraction is setting new benchmarks

More than 35 years after silica-based DNA and RNA isolation was first scientifically documented,¹ Analytik Jena is launching a global innovation in nucleic acid extraction. SmartExtraction significantly accelerates and considerably simplifies the entire procedure. Most notably, the technology accommodates the trend towards maximum process automation.

In order to provide users with maximum freedom when selecting materials, SmartExtraction was designed to be platform-independent. The technology can be used with all of Analytik Jena's pipetting systems, including InnuPure® C16 & C96, GeneTheatre, and CyBio® FeliX, and is simple to adapt for use with any liquid handling system.² The required laboratory equipment is reduced to a thermal shaker and a magnetrack for manual applications.

In addition to simplifying procedures, SmartExtraction is also superior to other technologies in terms of yield,

DNA quality, and efficiency criteria: Thanks to high binding capacities, large amounts of high-molecular DNA can be extracted with the appropriate starting materials.

Compared with magnetic particle technology used in conjunction with automated pipetting extraction systems, the new technology significantly increases the amount of extracted nucleic acids in many applications, while substantially reducing the processing time required.

That's not optimization – that's a quantum leap!



¹ Bert Vogelstein, David Gillespie; „Preparative and analytical purification of DNA from agarose“ Proc. Natl. Acad. Sci. USA; Vol. 76, No. 2, page 615619, February 1979; Biochemistry

² Pipetting systems with 1 mL pipetting heads

SmartExtraction

We Change the Way to Prep

No phenol/chloroform

No ion exchanger

No silica material and/or spin-filter columns

No silica or magnetic particle suspensions

SmartExtraction Combines the Best of Both Worlds

Analytik Jena's patented extraction chemistry (DC-Technology®) with an intelligent, modified surface material for the binding of nucleic acids. Other solid phases such as filter materials, magnetic nano- or micro-particles, and silica with all of their disadvantages are superfluous.

DC-Technology® meets Smart Surfaces

- No phenol/chloroform
- No ion exchanger
- No silica material and/or spin-filter columns
- No silica or magnetic particle suspensions


The core of the extraction chemistry (DC-Technology®) - as basis for SmartExtraction - is the combination of chaotropic and non-chaotropic salts with low ionic strengths, which enables the development of efficient lysis and binding buffers.

The Smart Modified Surface selectively binds the desired nucleic acids and then elutes them with high efficiency.

We're revolutionizing the process of nucleic acid extraction!

Your Benefits

- Exceptional yield, quality, and quantity
- Extraction of high-molecular weight DNA
- No centrifugation:
optimal handling
- No carryover of magnetic particles:
eluted nucleic acids of maximum quality
- No silica materials:
optimized binding capacity for nucleic acids



Minimal Effort for Maximum Efficiency

The entire purification process is limited to simple pipetting and incubation steps. There's no need to bind the DNA to magnetic particles or silica materials, nor is there any need for subsequent handling.

Automated SmartExtraction – concentrating on what is important

SmartExtraction has revolutionized the process of isolating and purifying nucleic acids, for instance, by making a number of steps, additional materials, and special equipment superfluous.

Developed by experts, the Smart Modified Surface ensures that the entire extraction process occurs in the pipette tip.

Using existing workflows as an example, that means that none of the steps, from lysis to DNA desorption, require special materials or equipment such as magnets, vacuum stations, or centrifuges. The entire analysis process is limited to simply pipetting up and down. There's no need to bind the DNA to magnetic particles or silica materials, thus eliminating time-consuming steps such as resuspension, separation, and centrifugation.

Manual SmartExtraction – the solution for small batches

The unique SmartExtraction technology does not just achieve excellent results on automated platforms. It simplifies the manual extraction of genomic DNA from various starting materials – requiring much less equipment than conventional solutions: Instead of a centrifuge, only a thermal shaker is needed.

Analytik Jena offers ready-to-use kits, including SmartExtraction tubes and granulates with Smart Modified Surfaces for a variety of application areas.

Your Benefits


- Special equipment and materials are superfluous
- Simple procedure, exceptionally easy to automate
- Can be used in liquid handling systems of all kinds²
- Optimized, ready-to-use kits for InnuPure[®] extraction systems and the CyBio[®] FeliX, and GeneTheatre liquid handling systems
- "Smart Modified Surface" for simple extraction in the pipette tip

Your Benefits

- Reduction of the required laboratory equipment
- Nucleic acid extraction in a thermal shaker
- Granulates with Smart Modified Surfaces in a SmartExtraction tube
- Quick, simple routines for the manual isolation of nucleic acids

The Results Speak for Themselves

In addition to a range of different starting materials, SmartExtraction is also suitable for samples with a high initial volume or high application amount in particular. All automated protocols are prepared in such a way that they don't require vacuum stations or magnetic adapters to work directly on pipetting platforms with a 1 ml filter tip.



Lane 1 and 14:

DNA ladder

Lane 2 – 5:

DNA extraction from the nucleated cells of 1 ml whole blood

Lane 6 – 9:

DNA extraction from the nucleated cells of 2 ml whole blood

Lane 10 – 13:

DNA extraction from the nucleated cells of 3 ml whole blood

The smart Blood DNA Midi prep (a) was used in fourfold determination for preparing 1 ml, 2 ml, and 3 ml whole blood samples. Following red blood cell lysis, the nucleated cells were pelleted and then resuspended in phosphate-buffered saline (PBS). All additional steps of DNA extraction were automated with InnuPure® C16.


Results table

Sample	Amount of whole blood [ml]*	Conc. [ng/μl]	Yield [μg]	$A_{260}:A_{280}$	$A_{260}:A_{230}$
1	1 mL	79	23.7	1.9	2.2
2	1 mL	71	21.3	1.9	2.2
3	1 mL	68	20.4	1.9	2.2
4	1 mL	63	18.9	1.9	2.2
5	2 mL	150	45.0	1.9	2.2
6	2 mL	156	46.8	1.9	2.2
7	2 mL	143	42.9	1.9	2.2
8	2 mL	147	44.1	1.9	2.3
9	3 mL	216	64.8	1.8	2.3
10	3 mL	234	70.2	1.9	2.3
11	3 mL	228	68.4	1.9	2.3
12	3 mL	240	72.0	1.9	2.3

* The resulting nucleated cells

With a comparable purity, the yield of the nucleic acid almost increases linearly with the amount of nucleated blood cells used. The gel image also clearly shows the exceptional quality of the high-molecular DNA extracted.

A comparable experiment isolated the DNA of gram-positive and negative bacteria. Using the SmartExtraction technology, extraction was carried out on five different liquid handling systems at the same time.



Gram-negative - *E. coli*


Device	Conc. [ng/μl]	Yield [μg]
InnuPure® C16	68.7	20.6
GeneTheatre	65.4	19.6
CyBio® FeliX	65.1	19.5
InnuPure® C96	65.7	19.7

Gram-positive - *S. aureus*

Device	Conc. [ng/μl]	Yield [μg]
InnuPure® C16	59.5	17.8
GeneTheatre	54.5	16.3
CyBio® FeliX	52.4	15.7
InnuPure® C96	49.5	14.8

The device comparison clearly reveals the universality of the extraction technology. High-quality nucleic acids are extracted with outstanding yield regardless of the device.

By the SmartExtraction tubes, the genomic DNA from different strains of gram- and gram+ bacteria were isolated in highly-efficient, manual processes.



Lane 1 and 10:

DNA ladders

Lane 2 to 9:

DNA extracted from different bacterial strains using SmartExtraction (each in double determination)

Results table

Lane	Description	Conc. [ng/μl]	$A_{260}:A_{280}$	$A_{260}:A_{230}$	Yield [μg]
1	Klebsiella pneumoniae	99.5	1.951	1.579	29.9
2	DSM 25721 (gram-)	102.0	2.040	2.147	30.6
3	Escherichia coli DSM 498	685.0	1.983	1.966	205.5
4	(gram-)	686.0	1.984	2.052	205.8
5	Enterococcus faecalis	240.0	1.996	1.886	72.0
6	DSM 20478 (gram+)	338.0	1.985	1.968	101.4
7	Staphylococcus aureus	144.0	1.901	1.597	43.2
8	BK916 (gram+)	134.0	1.942	1.914	40.2

The manual application of SmartExtraction technology also enables a very high yield of macromolecular DNA from both gram- and gram+ bacterial strains. The quality of the extracted DNA is excellent with regard to $A_{260}:A_{280}$ and $A_{260}:A_{230}$ assessment.

All on One View

SmartExtraction includes a number of ready-to-use kits for automated systems (a) as well as for manual (m) preparation of different starting materials.

Sample type	smart Blood DNA Midi prep (a)	smart DNA prep (a)	smart Blood DNA Midi Direct prep (a)	smart DNA prep (m)	smart Blood DNA Midi prep (m)
Starting material					
Whole Blood					
Up to 1 ml	✓		✓		✓
Up to 3 ml	✓				✓
Up to 10 ml					✓
Tissue					
Rodent tail pieces					
Mouse tail pieces		✓		✓	
Rat tail pieces		✓		✓	
Eucaryotic cells					
Bacteria (gram+ and gram-)		✓		✓	
Yeast cells		✓		✓	

Order Information

smart Blood DNA Midi prep (a)

Order number	Description	Compatible Analytik Jena devices
845-ASS-1208016	16 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASS-1208096	96 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-1208016	16 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-1208096	96 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-1296096	1 x 96 Reactions, Reagent Plates	InnuPure® C96, CyBio® FeliX

smart DNA prep (a)

Order number	Description	Compatible Analytik Jena devices
845-ASS-2008016	16 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASS-2008096	96 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-2008016	16 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-2008096	96 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, GeneTheatre, CyBio® FeliX
845-ASP-2096096	1 x 96 Reactions, Reagent Plates	InnuPure® C96, CyBio® FeliX

smart Blood DNA Midi Direct prep (a)

Order number	Description	Compatible Analytik Jena devices
845-ASS-3008016	16 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, CyBio® FeliX
845-ASS-3008096	96 Reactions, Reagent Strips, 1 Reaction per Strip	InnuPure® C16, CyBio® FeliX
845-ASP-3008016	16 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, CyBio® FeliX
845-ASP-3008096	96 Reactions, Reagent Plates, 8 Reactions per Plate	InnuPure® C16, CyBio® FeliX
845-ASP-3096096	1 x 96 Reactions, Reagent Plates	InnuPure® C96, CyBio® FeliX


smart DNA prep (m)

Order number	Description
845-KS-8000010	smart DNA prep (m), 10 Reactions
845-KS-8000050	smart DNA prep (m), 50 Reactions

smart Blood DNA Midi prep (m)

Order number	Description
845-KS-8100010	smart Blood DNA Midi prep (m), 10 Reactions
845-KS-8100050	smart Blood DNA Midi prep (m), 50 Reactions

You will find more information about our kits on our website www.analytik-jena.com.



Headquarters

Analytik Jena AG
Konrad-Zuse-Str. 1
07745 Jena · Germany

Phone +49 3641 77 70
Fax +49 3641 77 9279
info@analytik-jena.com
www.analytik-jena.com

Pictures: Analytik Jena AG
Subject to changes in design and scope of delivery as well as further technical development!